

# Alto™

Abdominal Stent  
Graft System


## STENT GRAFT SYSTEM OVERVIEW

### STEP 1

#### Select ALTO Aortic Body

Stent Graft Diameter (mm)	Aortic ID* (mm)
20	16-17
23	18-20
26	21-23
29	24-26
34	27-30


\*At the intended proximal sealing ring location.



### STEP 2


#### Select Ovation iX Iliac Limbs

Labeled Diameter (mm)	Native Iliac Vessel ID Range (mm)
10	8-9
12	10-11
14	12-13
16	14-15
18	16-17
22	18-20
28	21-25



### STEP 3

#### Select Ovation iX Extensions



## STENT GRAFT SYSTEM SPECIFICATIONS

ALTO AORTIC BODY	Catalog Number	Delivery System Inner Profile (F)	Delivery System Outer Profile (F)	Stent Graft Proximal Diameter (mm)	Stent Graft Length (mm)
	TV-AB2080-L	13	15	20	80
	TV-AB2380-L	13	15	23	80
	TV-AB2680-L	13	15	26	80
	TV-AB2980-L	13	15	29	80
	TV-AB3480-L	13	15	34	80

OVATION IX ILIAC LIMB	Catalog Number	Delivery System Inner Profile (F)	Delivery System Outer Profile (F)	Stent Graft Proximal Diameter (mm)	Stent Graft Distal Diameter (mm)	Stent Graft Length (mm)	Stent Graft Length with Aortic Body (mm)
	TV-IL141080-I	10	12	14	10	80	125
	TV-IL1410100-I	10	12	14	10	100	145
	TV-IL1410120-I	10	12	14	10	120	165
	TV-IL1410140-I	10	12	14	10	140	185
	TV-IL1410160-I	10	12	14	10	160	205
	TV-IL141280-I	10	12	14	12	80	125
	TV-IL1412100-I	10	12	14	12	100	145
	TV-IL1412120-I	10	12	14	12	120	165
	TV-IL1412140-I	10	12	14	12	140	185
	TV-IL1412160-I	10	12	14	12	160	205
	TV-IL141480-I	10	12	14	14	80	125
	TV-IL1414100-I	10	12	14	14	100	145
	TV-IL1414120-I	10	12	14	14	120	165
	TV-IL1414140-I	10	12	14	14	140	185
	TV-IL1414160-I	10	12	14	14	160	205
	TV-IL141680-I	11	13	14	16	80	125
	TV-IL1416100-I	11	13	14	16	100	145
	TV-IL1416120-I	11	13	14	16	120	165
	TV-IL1416140-I	11	13	14	16	140	185
	TV-IL1416160-I	11	13	14	16	160	205
	TV-IL141880-I	11	13	14	18	80	125
	TV-IL1418100-I	11	13	14	18	100	145
	TV-IL1418120-I	11	13	14	18	120	165
	TV-IL1418140-I	11	13	14	18	140	185
	TV-IL1418160-I	11	13	14	18	160	205
	TV-IL142280-I	12	14	14	22	80	125
	TV-IL1422100-I	12	14	14	22	100	145
	TV-IL1422120-I	12	14	14	22	120	165
	TV-IL1422140-I	12	14	14	22	140	185
	TV-IL1422160-I	12	14	14	22	160	205
	TV-IL142880-I	13	15	14	28	80	125
	TV-IL1428100-I	13	15	14	28	100	145
	TV-IL1428120-I	13	15	14	28	120	165
	TV-IL1428140-I	13	15	14	28	140	185
	TV-IL1428160-I	13	15	14	28	160	205

OVATION IX ILIAC EXTENSION	Catalog Number	Delivery System Inner Profile (F)	Delivery System Outer Profile (F)	Stent Graft Proximal Diameter (mm)	Stent Graft Distal Diameter (mm)	Stent Graft Length (mm)
	TV-EX101045-I	10	12	10	10	45
	TV-EX121245-I	10	12	12	12	45
	TV-EX141445-I	10	12	14	14	45
	TV-EX161645-I	11	13	16	16	45
	TV-EX181845-I	11	13	18	18	45
	TV-EX222245-I	12	14	22	22	45
	TV-EX282845-I	13	15	28	28	45

CUSTOMSEAL™ POLYMER FILL KIT	Catalog Number
	TV-CS14-F

AUTOINJECTOR 2	Catalog Number
	TV-AI01-L



The ALTO™ Abdominal Stent Graft System and associated components are not available in all countries or regions.

Please contact your Endologix representative for details regarding product availability. Prior to use, refer to Instructions for Use for more information concerning Indications, Contraindications, Specific Anatomic Considerations, Warnings, Precautions, and Adverse Events.

Rx only.

©2020 Endologix Inc. All rights reserved. MM2146 Rev 03

Not for distribution in the USA.

