


Endurant™ II/IIIs

AAA stent graft system

Endurant™ II/IIIs system product code description

ET	B	F	23	13	C	124	E	18
								Catheter outer diameter
							Delivery system E - Endurant™ II system	
						Total covered length		
					Distal design C - Closed web			
				Distal graft diameter				
			Proximal graft diameter					
		Proximal design F - FreeFlo W - Open web						
		Device configuration B - Bifurcations L - Limbs E - Iliac extension C - Extensions & cuffs T - Extensions & cuffs U - Aorto-uni-iliac (AUI)						
		Product name ET - Endurant™ II system ES - Endurant™ IIs system						



Endurant™ IIs system bifurcations

	Product code					
	Proximal graft diameter (mm)	Distal graft diameter (mm)	Distal design	Total covered length (mm)	Delivery system	Catheter outer diameter (Fr)
ESBF	23	14	C	103	E	18
ESBF	25	14	C	103	E	18
ESBF	28	14	C	103	E	18
ESBF	32	14	C	103	E	20
ESBF	36	14	C	103	E	20



Endurant™ II system bifurcations

Product code						
	Proximal graft diameter (mm)	Distal graft diameter (mm)	Distal design	Total covered length (mm)	Delivery system	Catheter outer diameter (Fr)
ETBF	23	13	C	124	E	18
ETBF	23	13	C	145	E	18
ETBF	23	13	C	166	E	18
ETBF	23	16	C	124	E	18
ETBF	23	16	C	145	E	18
ETBF	23	16	C	166	E	18
ETBF	25	13	C	124	E	18
ETBF	25	13	C	145	E	18
ETBF	25	13	C	166	E	18
ETBF	25	16	C	124	E	18
ETBF	25	16	C	145	E	18
ETBF	25	16	C	166	E	18
ETBF	28	13	C	124	E	18
ETBF	28	13	C	145	E	18
ETBF	28	13	C	166	E	18
ETBF	28	16	C	124	E	18
ETBF	28	16	C	145	E	18
ETBF	28	16	C	166	E	18
ETBF	28	20	C	124	E	18
ETBF	28	20	C	145	E	18
ETBF	28	20	C	166	E	18
ETBF	32	16	C	124	E	20
ETBF	32	16	C	145	E	20
ETBF	32	16	C	166	E	20
ETBF	32	20	C	124	E	20
ETBF	32	20	C	145	E	20
ETBF	32	20	C	166	E	20
ETBF	36	16	C	145	E	20
ETBF	36	16	C	166	E	20
ETBF	36	20	C	145	E	20
ETBF	36	20	C	166	E	20

Endurant™ II/IIs

AAA stent graft system

Limbs†

	Product code					Catheter outer diameter (Fr)	Total contralateral covered length with EII/EIIs bifurcated**	Total ipsilateral covered length with EIs bifurcated††
	Proximal graft diameter (mm)	Distal graft diameter (mm)	Distal design	Total covered length (mm)	Delivery system			
ETLW	16	10	C	82	E	14	132	152
ETLW	16	10	C	93	E	14	143	163
ETLW	16	10	C	124	E	14	174	174 - 194
ETLW	16	10	C	146	E	16	196	196 - 216
ETLW	16	10	C	156	E	16	206	206 - 226
ETLW	16	10	C	199	E	16	249	249 - 269
ETLW	16	13	C	82	E	14	132	152
ETLW	16	13	C	93	E	14	143	163
ETLW	16	13	C	124	E	14	174	174 - 194
ETLW	16	13	C	146	E	16	196	196 - 216
ETLW	16	13	C	156	E	16	206	206 - 226
ETLW	16	13	C	199	E	16	249	249 - 269
ETLW	16	16	C	82	E	14	132	132 - 152
ETLW	16	16	C	93	E	14	143	143 - 163
ETLW	16	16	C	124	E	14	174	174 - 194
ETLW	16	16	C	146	E	16	196	196 - 216
ETLW	16	16	C	156	E	16	206	206 - 226
ETLW	16	16	C	199	E	16	249	249 - 269
ETLW	16	20	C	82	E	16	132	152
ETLW	16	20	C	93	E	16	143	163
ETLW	16	20	C	124	E	16	174	174 - 194
ETLW	16	20	C	146	E	16	196	196 - 216
ETLW	16	20	C	156	E	16	206	206 - 226
ETLW	16	20	C	199	E	16	249	249 - 269
ETLW	16	24	C	82	E	16	132	152
ETLW	16	24	C	93	E	16	143	163
ETLW	16	24	C	124	E	16	174	174 - 194
ETLW	16	24	C	146	E	16	196	196 - 216
ETLW	16	24	C	156	E	16	206	206 - 226
ETLW	16	24	C	199	E	16	249	249 - 269
ETLW	16	28	C	82	E	16	132	152
ETLW	16	28	C	93	E	16	143	163
ETLW	16	28	C	124	E	16	174	174 - 194
ETLW	16	28	C	146	E	16	196	196 - 216
ETLW	16	28	C	156	E	16	206	206 - 226
ETLW	16	28	C	199	E	16	249	249 - 269

† The limb mates with the AUI stent graft on the ipsilateral side.

** These calculations assume the minimum 30 mm overlap between the bifurcated stent graft and the contralateral iliac limb per the Endurant™ II stent graft system *Instructions for Use*. When using the 124 mm length bifurcated stent graft, subtract 10 mm from total contralateral covered length with bifurcated.

†† The 3-5 stent overlap is available only with select limbs. Please refer to the *Instructions for Use* for more information.



Iliac extensions

Product code						
	Proximal graft diameter (mm)	Distal graft diameter (mm)	Distal design	Total covered length (mm)	Delivery system	Catheter outer diameter (Fr)
ETEW	10	10	C	82	E	14
ETEW	13	13	C	82	E	14
ETEW	20	20	C	82	E	16
ETEW	24	24	C	82	E	16
ETEW	28	28	C	82	E	18

Aortic extensions

Product code						
	Proximal graft diameter (mm)	Distal graft diameter (mm)	Distal design	Total covered length (mm)	Delivery system	Catheter outer diameter (Fr)
ETCF	23	23	C	49	E	18
ETCF	25	25	C	49	E	18
ETCF	28	28	C	49	E	18
ETCF	32	32	C	49	E	20
ETCF	36	36	C	49	E	20
ETTF	23	23	C	70	E	18
ETTF	25	25	C	70	E	18
ETTF	28	28	C	70	E	18
ETTF	32	32	C	70	E	20
ETTF	36	36	C	70	E	20


AUI

Product code						
	Proximal graft diameter (mm)	Distal graft diameter (mm)	Distal design	Total covered length (mm)	Delivery system	Catheter outer diameter (Fr)
ETUF	23	14	C	102	E	18
ETUF	25	14	C	102	E	18
ETUF	28	14	C	102	E	18
ETUF	32	14	C	102	E	20
ETUF	36	14	C	102	E	20



Endurant™ II system bifurcations

	Product code					Graft dimensions (mm)					
	Proximal graft diameter (mm)	Distal graft diameter (mm)	Distal design	Total covered length (mm)	Delivery system	A	B	C	D	E	F
ETBF	23	13	C	124	E	12	40	74	25	-	-
ETBF	23	13	C	145	E	12	50	84	25	-	-
ETBF	23	13	C	166	E	12	50	84	25	-	-
ETBF	23	16	C	124	E	12	40	74	25	10	30
ETBF	23	16	C	145	E	12	50	84	25	10	40
ETBF	23	16	C	166	E	12	50	84	25	10	60
ETBF	25	13	C	124	E	14	40	74	27	-	-
ETBF	25	13	C	145	E	14	50	84	27	-	-
ETBF	25	13	C	166	E	14	50	84	27	-	-
ETBF	25	16	C	124	E	14	40	74	30	-	-
ETBF	25	16	C	145	E	14	50	84	30	-	-
ETBF	25	16	C	166	E	14	50	84	30	-	-
ETBF	28	13	C	124	E	14	40	74	27	-	-
ETBF	28	13	C	145	E	14	50	84	27	-	-
ETBF	28	13	C	166	E	14	50	84	27	-	-
ETBF	28	16	C	124	E	14	40	74	30	-	-
ETBF	28	16	C	145	E	14	50	84	30	-	-
ETBF	28	16	C	166	E	14	50	84	30	-	-
ETBF	28	20	C	124	E	14	40	74	30	10	30
ETBF	28	20	C	145	E	14	50	84	30	10	40
ETBF	28	20	C	166	E	14	50	84	30	10	60
ETBF	32	16	C	124	E	14	40	74	30	-	-
ETBF	32	16	C	145	E	14	50	84	30	-	-
ETBF	32	16	C	166	E	14	50	84	30	-	-
ETBF	32	20	C	124	E	14	40	74	30	10	30
ETBF	32	20	C	145	E	14	50	84	30	10	40
ETBF	32	20	C	166	E	14	50	84	30	10	60
ETBF	36	16	C	145	E	14	50	84	30	-	-
ETBF	36	16	C	166	E	14	50	84	30	-	-
ETBF	36	20	C	145	E	14	50	84	30	10	40
ETBF	36	20	C	166	E	14	50	84	30	10	60



Endurant™ IIs system bifurcations


Product code						Graft dimensions (mm)			
	Proximal graft diameter (mm)	Distal graft diameter (mm)	Distal design	Total covered length (mm)	Delivery system	A	B	C	D
ESBF	23	14	C	103	E	14	50	84	28
ESBF	25	14	C	103	E	14	50	84	28
ESBF	28	14	C	103	E	14	50	84	28
ESBF	32	14	C	103	E	14	50	84	28
ESBF	36	14	C	103	E	14	50	84	28

Endurant™ II/IIIs

AAA stent graft system

Tapered limbs

	Product code					Graft dimensions (mm)	
	Proximal graft diameter (mm)	Distal graft diameter (mm)	Distal design	Total covered length (mm)	Delivery system	A	B
ETLW	16	10	C	82	E	20	30
ETLW	16	10	C	93	E	20	40
ETLW	16	10	C	124	E	20	40
ETLW	16	10	C	146	E	20	62
ETLW	16	10	C	156	E	20	72
ETLW	16	10	C	199	E	20	115
ETLW	16	13	C	82	E	10	40
ETLW	16	13	C	93	E	10	40
ETLW	16	13	C	124	E	10	40
ETLW	16	13	C	146	E	10	62
ETLW	16	13	C	156	E	10	72
ETLW	16	13	C	199	E	10	115




Tapered limbs



Straight limbs

Product code					
	Proximal graft diameter (mm)	Distal graft diameter (mm)	Distal design	Total covered length (mm)	Delivery system
ETLW	16	16	C	82	E
ETLW	16	16	C	93	E
ETLW	16	16	C	124	E
ETLW	16	16	C	146	E
ETLW	16	16	C	156	E
ETLW	16	16	C	199	E



Straight limbs



Flared limbs


	Product code					Graft dimensions (mm)	
	Proximal graft diameter (mm)	Distal graft diameter (mm)	Distal design	Total covered length (mm)	Delivery system	A	B
	ETLW	16	20	C	82	E	10
ETLW	16	20	C	93	E	10	40
ETLW	16	20	C	124	E	10	40
ETLW	16	20	C	146	E	10	40
ETLW	16	20	C	156	E	10	40
ETLW	16	20	C	199	E	10	40
ETLW	16	24	C	82	E	20	30
ETLW	16	24	C	93	E	20	40
ETLW	16	24	C	124	E	20	40
ETLW	16	24	C	146	E	20	40
ETLW	16	24	C	156	E	20	40
ETLW	16	24	C	199	E	20	40
ETLW	16	28	C	82	E	20	30
ETLW	16	28	C	93	E	20	40
ETLW	16	28	C	124	E	20	40
ETLW	16	28	C	146	E	20	40
ETLW	16	28	C	156	E	20	40
ETLW	16	28	C	199	E	20	40

Endurant™ II/IIIs

AAA stent graft system

Iliac extensions

Product code					
	Proximal graft diameter (mm)	Distal graft diameter (mm)	Distal design	Total covered length (mm)	Delivery system
ETEW	10	10	C	82	E
ETEW	13	13	C	82	E
ETEW	20	20	C	82	E
ETEW	24	24	C	82	E
ETEW	28	28	C	82	E




Iliac extensions



AUI


	Product code						Graft dimensions (mm)	
	Proximal graft diameter (mm)	Distal graft diameter (mm)	Distal design	Total covered length (mm)	Delivery system	A	B	
ETUF	23	14	C	102	E	30	40	
ETUF	25	14	C	102	E	30	40	
ETUF	28	14	C	102	E	30	40	
ETUF	32	14	C	102	E	30	40	
ETUF	36	14	C	102	E	30	40	

The only device with an FDA-approved AUI Indication†



Endurant™ II AUI

† As of November 2020



Aortic extensions

	Product code				
	Proximal graft diameter (mm)	Distal graft diameter (mm)	Distal design	Total covered length (mm)	Delivery system
ETCF	23	23	C	49	E
ETCF	25	25	C	49	E
ETCF	28	28	C	49	E
ETCF	32	32	C	49	E
ETCF	36	36	C	49	E
ETTF	23	23	C	70	E
ETTF	25	25	C	70	E
ETTF	28	28	C	70	E
ETTF	32	32	C	70	E
ETTF	36	36	C	70	E