

Time-tested
accuracy and
durability.

Zenith Flex[®]

AAA ENDOVASCULAR GRAFT

Bifurcated Main Body Grafts

Order Number	Reference Part Number	Aortic Vessel OD mm	Main Body Diameter mm	Main Body Contralateral Length mm	Introducer Sheath ID/OD Fr (mm)/mm
22 mm Diameter					
G48402	TFFB-22-82-ZT	18-19	22	82	18 (6.0)/7.1
G48403	TFFB-22-96-ZT	18-19	22	96	18 (6.0)/7.1
G48404	TFFB-22-111-ZT	18-19	22	111	18 (6.0)/7.1
G48405	TFFB-22-125-ZT	18-19	22	125	18 (6.0)/7.1
G48406	TFFB-22-140-ZT	18-19	22	140	18 (6.0)/7.1
24 mm Diameter					
G48407	TFFB-24-82-ZT	20-21	24	82	18 (6.0)/7.1
G48408	TFFB-24-96-ZT	20-21	24	96	18 (6.0)/7.1
G48409	TFFB-24-111-ZT	20-21	24	111	18 (6.0)/7.1
G48410	TFFB-24-125-ZT	20-21	24	125	18 (6.0)/7.1
G48411	TFFB-24-140-ZT	20-21	24	140	18 (6.0)/7.1
26 mm Diameter					
G48412	TFFB-26-82-ZT	22	26	82	18 (6.0)/7.1
G48413	TFFB-26-96-ZT	22	26	96	18 (6.0)/7.1
G48414	TFFB-26-111-ZT	22	26	111	18 (6.0)/7.1
G48415	TFFB-26-125-ZT	22	26	125	18 (6.0)/7.1
G48416	TFFB-26-140-ZT	22	26	140	18 (6.0)/7.1
28 mm Diameter					
G48417	TFFB-28-82-ZT	23-24	28	82	20 (6.7)/7.7
G48418	TFFB-28-96-ZT	23-24	28	96	20 (6.7)/7.7
G48419	TFFB-28-111-ZT	23-24	28	111	20 (6.7)/7.7
G48420	TFFB-28-125-ZT	23-24	28	125	20 (6.7)/7.7
G48421	TFFB-28-140-ZT	23-24	28	140	20 (6.7)/7.7
30 mm Diameter					
G48422	TFFB-30-82-ZT	25-26	30	82	20 (6.7)/7.7
G48423	TFFB-30-96-ZT	25-26	30	96	20 (6.7)/7.7
G48424	TFFB-30-111-ZT	25-26	30	111	20 (6.7)/7.7
G48425	TFFB-30-125-ZT	25-26	30	125	20 (6.7)/7.7
G48426	TFFB-30-140-ZT	25-26	30	140	20 (6.7)/7.7

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Bifurcated Main Body Grafts (Continued)

Order Number	Reference Part Number	Aortic Vessel OD mm	Main Body Diameter mm	Main Body Contralateral Length mm	Introducer Sheath ID/OD Fr (mm)/mm
32 mm Diameter					
G48427	TFFB-32-82-ZT	27-28	32	82	20 (6.7)/7.7
G48428	TFFB-32-96-ZT	27-28	32	96	20 (6.7)/7.7
G48429	TFFB-32-111-ZT	27-28	32	111	20 (6.7)/7.7
G48430	TFFB-32-125-ZT	27-28	32	125	20 (6.7)/7.7
G48431	TFFB-32-140-ZT	27-28	32	140	20 (6.7)/7.7
36 mm Diameter					
G48432	TFFB-36-95-ZT	29-32	36	95	22 (7.3)/8.5
G48433	TFFB-36-113-ZT	29-32	36	113	22 (7.3)/8.5
G48434	TFFB-36-131-ZT	29-32	36	131	22 (7.3)/8.5
G48435	TFFB-36-149-ZT	29-32	36	149	22 (7.3)/8.5


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Z-Trak® Introduction System for Main Body Grafts

22-26 mm grafts: 18 Fr (6.0 mm) ID/7.1 mm OD

28-32 mm grafts: 20 Fr (6.7 mm) ID/7.7 mm OD

36 mm grafts: 22 Fr (7.3 mm) ID/8.5 mm OD



Ancillary Components

Order Number	Reference Part Number	Diameter mm	Length mm	Introducer Sheath ID/OD Fr (mm)/mm
Main Body Extensions				
G55204	ESBE-22-39-ZT	22	39	18 (6.0)/7.1
G55205	ESBE-24-39-ZT	24	39	18 (6.0)/7.1
G55206	ESBE-26-39-ZT	26	39	18 (6.0)/7.1
G55207	ESBE-28-39-ZT	28	39	20 (6.7)/7.7
G55208	ESBE-30-39-ZT	30	39	20 (6.7)/7.7
G55209	ESBE-32-39-ZT	32	39	20 (6.7)/7.7
G55210	ESBE-36-50-ZT	36	50	20 (6.7)/7.7

Order Number	Reference Part Number	Diameter mm	Length mm	Introducer Sheath ID/OD Fr (mm)/mm
Main Body Extensions				
G55211	ESBE-22-58-ZT	22	58	18 (6.0)/7.1
G55212	ESBE-24-58-ZT	24	58	18 (6.0)/7.1
G55213	ESBE-26-58-ZT	26	58	18 (6.0)/7.1
G55214	ESBE-28-58-ZT	28	58	20 (6.7)/7.7
G55215	ESBE-30-58-ZT	30	58	20 (6.7)/7.7
G55216	ESBE-32-58-ZT	32	58	20 (6.7)/7.7
G55217	ESBE-36-73-ZT	36	73	20 (6.7)/7.7

Order Number	Reference Part Number	Proximal Diameter mm	Distal Diameter mm	Length mm	Introducer Sheath ID/OD Fr (mm)/mm
Converters					
G55191	ESC-24-12-80-ZT	24	12	80	18 (6.0)/7.1
G55192	ESC-28-12-80-ZT	28	12	80	20 (6.7)/7.7
G55193	ESC-32-12-80-ZT	32	12	80	20 (6.7)/7.7
G55194	ESC-36-12-82-ZT	36	12	82	20 (6.7)/7.7


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Ancillary Components (Continued)

Order Number	Reference Part Number	Diameter mm	Length mm	Introducer Sheath ID/OD Fr (mm)/mm
Iliac Plugs				
G18505	ZIP-14-30-ZT	14	30	14 (4.7)/5.4
G18506	ZIP-16-30-ZT	16	30	14 (4.7)/5.4
G18507	ZIP-20-30-ZT	20	30	16 (5.3)/6.0
G18508	ZIP-24-30-ZT	24	30	16 (5.3)/6.0

Order Number	Reference Part Number	Proximal Diameter mm	Distal Diameter mm	Length mm	Introducer Sheath ID/OD Fr (mm)/mm
Iliac Leg Extension					
G55197	ESLE-12-55-ZT	12	12	55	14 (4.7)/5.4
G55198	ESLE-14-55-ZT	14	14	55	14 (4.7)/5.4
G55199	ESLE-16-55-ZT	16	16	55	14 (4.7)/5.4
G55200	ESLE-18-55-ZT	18	18	55	16 (5.3)/6.0
G55201	ESLE-20-55-ZT	20	20	55	16 (5.3)/6.0
G55202	ESLE-22-55-ZT	22	22	55	16 (5.3)/6.0
G55203	ESLE-24-55-ZT	24	24	55	16 (5.3)/6.0

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.




Helps restore
graft-to-vessel
fixation and seal.

Zenith Renu[®]

AAA ANCILLARY GRAFT

Indications for Use



Indications for Use	Secondary endovascular intervention in patients having received prior endovascular repair of infrarenal abdominal aortic or aortoiliac aneurysms in which there is inadequate proximal fixation or seal.		
Length	Lowest renal artery to bifurcation of previously placed graft of > 43 mm for the main body extension and > 37 mm for the converter.		
Aortic Fixation Site Diameter	18-32 mm (measured outer wall to outer wall)		
Angulation/ Curvature	<ul style="list-style-type: none"> • Angle < 60 degrees relative to the long axis of the aneurysm • Angle < 45 degrees relative to the axis of the suprarenal aorta 		
Converter Distal Fixation Site	For the converter used without an iliac leg, distal fixation site within a graft segment of ≤ 12 mm in diameter and ≥ 17 mm (one Cook-Z [®] stent) in length, with more overlap length being preferred. For the converter used in combination with the iliac leg, distal fixation site 7.5-20 mm in diameter (measured outer wall to outer wall) and > 10 mm in length with 20 to 30 mm being preferred.		
Main Body Distal Fixation Site	Main body extension distal fixation site within a graft segment of ≤ 34 mm diameter and at least one Cook-Z stent in length, with more overlap length being preferred.		
Access	Adequate femoral/iliac access compatible with the required introduction systems		
Introducer Diameters	22-26 mm grafts: ID 18 Fr (6.0 mm) OD 7.1 mm	28-32 mm grafts: ID 20 Fr (6.7 mm) OD 7.7 mm	36 mm grafts: ID 22 Fr (7.3 mm) OD 8.5 mm
Contraindications	<ul style="list-style-type: none"> • Patients with known sensitivities or allergies to stainless steel, polyester, solder (tin, silver), nitinol, polypropylene or gold. • Patients with a systemic infection who may be at increased risk of endovascular graft infection. 		

Converters

Order Number	Reference Part Number	Aortic Vessel OD mm	Proximal/Distal Diameter mm	Proximal/Distal Length mm	Total Length mm	Introducer Sheath ID/OD Fr (mm)/mm
22 mm Diameter						
G55355	AX1-1-22-108-ZT	18-19	22/12	37/71	108	18 (6.0)/7.1
G55356	AX1-1-22-125-ZT	18-19	22/12	37/88	125	18 (6.0)/7.1
G55357	AX1-2-22-113-ZT	18-19	22/12	59/54	113	18 (6.0)/7.1
24 mm Diameter						
G55361	AX1-1-24-108-ZT	20-21	24/12	37/71	108	18 (6.0)/7.1
G55362	AX1-1-24-125-ZT	20-21	24/12	37/88	125	18 (6.0)/7.1
G55363	AX1-2-24-113-ZT	20-21	24/12	59/54	113	18 (6.0)/7.1
26 mm Diameter						
G55367	AX1-1-26-108-ZT	22	26/12	37/71	108	18 (6.0)/7.1
G55368	AX1-1-26-125-ZT	22	26/12	37/88	125	18 (6.0)/7.1
G55369	AX1-2-26-113-ZT	22	26/12	59/54	113	18 (6.0)/7.1
28 mm Diameter						
G55374	AX1-1-28-108-ZT	23-24	28/12	37/71	108	20 (6.7)/7.7
G55375	AX1-1-28-125-ZT	23-24	28/12	37/88	125	20 (6.7)/7.7
G55376	AX1-2-28-113-ZT	23-24	28/12	59/54	113	20 (6.7)/7.7
30 mm Diameter						
G55380	AX1-1-30-108-ZT	25-26	30/12	37/71	108	20 (6.7)/7.7
G55381	AX1-1-30-125-ZT	25-26	30/12	37/88	125	20 (6.7)/7.7
G55382	AX1-2-30-113-ZT	25-26	30/12	59/54	113	20 (6.7)/7.7
32 mm Diameter						
G55386	AX1-1-32-108-ZT	27-28	32/12	37/71	108	20 (6.7)/7.7
G55387	AX1-1-32-125-ZT	27-28	32/12	37/88	125	20 (6.7)/7.7
G55388	AX1-2-32-113-ZT	27-28	32/12	59/54	113	20 (6.7)/7.7
36 mm Diameter						
G55392	AX1-1-36-116-ZT	29-32	36/12	52/64	116	22 (7.3)/8.5
G55394	AX1-2-36-127-ZT	29-32	36/12	80/47	127	22 (7.3)/8.5


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Z-Trak® Introduction System for Renu Converters

22-26 mm grafts: 18 Fr (6.0 mm) ID/7.1 mm OD

28-32 mm grafts: 20 Fr (6.7 mm) ID/7.7 mm OD

36 mm grafts: 22 Fr (7.3 mm) ID/8.5 mm OD



Main Body Extensions

Order Number	Reference Part Number	Aortic Vessel OD mm	Diameter mm	Length mm	Introducer Sheath ID/OD Fr (mm)/7.1
22 mm Diameter					
G55397	RX1-22-43-ZT	18-19	22	43	18 (6.0)/7.1
G55398	RX1-22-62-ZT	18-19	22	62	18 (6.0)/7.1
24 mm Diameter					
G55400	RX1-24-43-ZT	20-21	24	43	18 (6.0)/7.1
G55401	RX1-24-62-ZT	20-21	24	62	18 (6.0)/7.1
26 mm Diameter					
G55403	RX1-26-43-ZT	22	26	43	18 (6.0)/7.1
G55404	RX1-26-62-ZT	22	26	62	18 (6.0)/7.1
28 mm Diameter					
G55406	RX1-28-43-ZT	23-24	28	43	20 (6.7)/7.7
G55407	RX1-28-62-ZT	23-24	28	62	20 (6.7)/7.7
30 mm Diameter					
G55409	RX1-30-43-ZT	25-26	30	43	20 (6.7)/7.7
G55410	RX1-30-62-ZT	25-26	30	62	20 (6.7)/7.7
32 mm Diameter					
G55412	RX1-32-43-ZT	27-28	32	43	20 (6.7)/7.7
G55413	RX1-32-62-ZT	27-28	32	62	20 (6.7)/7.7
36 mm Diameter					
G55415	RX1-36-54-ZT	29-32	36	54	22 (7.3)/8.5
G55416	RX1-36-77-ZT	29-32	36	77	22 (7.3)/8.5


Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.


Z-Trak[®] Introduction System for Renu Main Body Extensions

22-26 mm grafts: 18 Fr (6.0 mm) ID/7.1 mm OD

28-32 mm grafts: 20 Fr (6.7 mm) ID/7.7 mm OD

36 mm grafts: 22 Fr (7.3 mm) ID/8.5 mm OD





Helps deliver
a durable repair
in complex
infrarenal anatomy.

Zenith[®] Fenestrated

AAA ENDOVASCULAR GRAFT


Proximal Body Grafts

Order Number	Reference Part Number	Aortic Vessel OD mm	Proximal Body Diameter mm	Proximal Body Length mm	Introducer Sheath ID/OD Fr (mm)/mm
One Proximal Internal Stent					
24 mm Diameter					
G32502	ZFEN-P-1-24-76-R	19-21	24	76	20 (6.7)/7.7
G32507	ZFEN-P-1-24-91-R	19-21	24	91	20 (6.7)/7.7
G32512	ZFEN-P-1-24-106-R	19-21	24	106	20 (6.7)/7.7
G32517	ZFEN-P-1-24-121-R	19-21	24	121	20 (6.7)/7.7
26 mm Diameter					
G32503	ZFEN-P-1-26-76-R	21-22	26	76	20 (6.7)/7.7
G32508	ZFEN-P-1-26-91-R	21-22	26	91	20 (6.7)/7.7
G32513	ZFEN-P-1-26-106-R	21-22	26	106	20 (6.7)/7.7
G32518	ZFEN-P-1-26-121-R	21-22	26	121	20 (6.7)/7.7
28 mm Diameter					
G32504	ZFEN-P-1-28-76-R	23-24	28	76	20 (6.7)/7.7
G32509	ZFEN-P-1-28-91-R	23-24	28	91	20 (6.7)/7.7
G32514	ZFEN-P-1-28-106-R	23-24	28	106	20 (6.7)/7.7
G32519	ZFEN-P-1-28-121-R	23-24	28	121	20 (6.7)/7.7
30 mm Diameter					
G32505	ZFEN-P-1-30-76-R	24-26	30	76	20 (6.7)/7.7
G32510	ZFEN-P-1-30-91-R	24-26	30	91	20 (6.7)/7.7
G32515	ZFEN-P-1-30-106-R	24-26	30	106	20 (6.7)/7.7
G32520	ZFEN-P-1-30-121-R	24-26	30	121	20 (6.7)/7.7
32 mm Diameter					
G32506	ZFEN-P-1-32-76-R	26-28	32	76	20 (6.7)/7.7
G32511	ZFEN-P-1-32-91-R	26-28	32	91	20 (6.7)/7.7
G32516	ZFEN-P-1-32-106-R	26-28	32	106	20 (6.7)/7.7
G32521	ZFEN-P-1-32-121-R	26-28	32	121	20 (6.7)/7.7
34 mm Diameter					
G32522	ZFEN-P-1-34-84-R	28-29	34	84	20 (6.7)/7.7
G32524	ZFEN-P-1-34-99-R	28-29	34	99	20 (6.7)/7.7
G32526	ZFEN-P-1-34-114-R	28-29	34	114	20 (6.7)/7.7
G32528	ZFEN-P-1-34-129-R	28-29	34	129	20 (6.7)/7.7
36 mm Diameter					
G32523	ZFEN-P-1-36-84-R	29-31	36	84	22 (7.3)/8.5
G32525	ZFEN-P-1-36-99-R	29-31	36	99	22 (7.3)/8.5
G32527	ZFEN-P-1-36-114-R	29-31	36	114	22 (7.3)/8.5
G32529	ZFEN-P-1-36-129-R	29-31	36	129	22 (7.3)/8.5

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

H&L-B One-Shot® Introduction System for Proximal Body Grafts

20 Fr (6.7 mm) ID/7.7 mm OD, 22 Fr (7.3 mm) ID/8.5 mm OD



Proximal Body Grafts


Order Number	Reference Part Number	Aortic Vessel OD mm	Proximal Body Diameter mm	Proximal Body Length mm	Introducer Sheath ID/OD Fr (mm)/mm
Two Proximal Internal Stents					
24 mm Diameter					
G32530	ZFEN-P-2-24-94-R	19-21	24	94	20 (6.7)/7.7
G32535	ZFEN-P-2-24-109-R	19-21	24	109	20 (6.7)/7.7
G32540	ZFEN-P-2-24-124-R	19-21	24	124	20 (6.7)/7.7
26 mm Diameter					
G32531	ZFEN-P-2-26-94-R	21-22	26	94	20 (6.7)/7.7
G32536	ZFEN-P-2-26-109-R	21-22	26	109	20 (6.7)/7.7
G32541	ZFEN-P-2-26-124-R	21-22	26	124	20 (6.7)/7.7
28 mm Diameter					
G32532	ZFEN-P-2-28-94-R	23-24	28	94	20 (6.7)/7.7
G32537	ZFEN-P-2-28-109-R	23-24	28	109	20 (6.7)/7.7
G32542	ZFEN-P-2-28-124-R	23-24	28	124	20 (6.7)/7.7
30 mm Diameter					
G32533	ZFEN-P-2-30-94-R	24-26	30	94	20 (6.7)/7.7
G32538	ZFEN-P-2-30-109-R	24-26	30	109	20 (6.7)/7.7
G32543	ZFEN-P-2-30-124-R	24-26	30	124	20 (6.7)/7.7
32 mm Diameter					
G32534	ZFEN-P-2-32-94-R	26-28	32	94	20 (6.7)/7.7
G32539	ZFEN-P-2-32-109-R	26-28	32	109	20 (6.7)/7.7
G32544	ZFEN-P-2-32-124-R	26-28	32	124	20 (6.7)/7.7
34 mm Diameter					
G32545	ZFEN-P-2-34-107-R	28-29	34	107	20 (6.7)/7.7
G32547	ZFEN-P-2-34-122-R	28-29	34	122	20 (6.7)/7.7
G32549	ZFEN-P-2-34-137-R	28-29	34	137	20 (6.7)/7.7
36 mm Diameter					
G32546	ZFEN-P-2-36-107-R	29-31	36	107	22 (7.3)/8.5
G32548	ZFEN-P-2-36-122-R	29-31	36	122	22 (7.3)/8.5
G32550	ZFEN-P-2-36-137-R	29-31	36	137	22 (7.3)/8.5

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

H&L-B One-Shot® Introduction System for Proximal Body Grafts

20 Fr (6.7 mm) ID/7.7 mm OD

22 Fr (7.3 mm) ID/8.5 mm OD



Distal Bifurcated Body Grafts

Order Number	Reference Part Number	Ipsilateral Leg Diameter mm	Distal Body Length mm	Distal Leg Length mm	Introducer Sheath ID/OD Fr (mm)/mm
Distal Bifurcated Body Graft					
G32551	ZFEN-D-12-28-76-C	12	76	28	20 (6.7)/7.7
G32554	ZFEN-D-12-28-91-C	12	91	28	20 (6.7)/7.7
G32500	ZFEN-D-12-28-106-C	12	106	28	20 (6.7)/7.7
G32501	ZFEN-D-12-28-121-C	12	121	28	20 (6.7)/7.7
G32561	ZFEN-D-12-45-76-C	12	76	45	20 (6.7)/7.7
G32565	ZFEN-D-12-45-91-C	12	91	45	20 (6.7)/7.7
G32569	ZFEN-D-12-45-106-C	12	106	45	20 (6.7)/7.7
G32573	ZFEN-D-12-45-121-C	12	121	45	20 (6.7)/7.7
G32577	ZFEN-D-12-62-76-C	12	76	62	20 (6.7)/7.7
G32580	ZFEN-D-12-62-91-C	12	91	62	20 (6.7)/7.7
G32584	ZFEN-D-12-62-106-C	12	106	62	20 (6.7)/7.7
G32588	ZFEN-D-12-62-121-C	12	121	62	20 (6.7)/7.7
G32552	ZFEN-D-16-28-76-C	16	76	28	20 (6.7)/7.7
G32555	ZFEN-D-16-28-91-C	16	91	28	20 (6.7)/7.7
G32557	ZFEN-D-16-28-106-C	16	106	28	20 (6.7)/7.7
G32559	ZFEN-D-16-28-121-C	16	121	28	20 (6.7)/7.7
G32562	ZFEN-D-16-45-76-C	16	76	45	20 (6.7)/7.7
G32566	ZFEN-D-16-45-91-C	16	91	45	20 (6.7)/7.7
G32570	ZFEN-D-16-45-106-C	16	106	45	20 (6.7)/7.7
G32574	ZFEN-D-16-45-121-C	16	121	45	20 (6.7)/7.7
G32592	ZFEN-D-16-62-76-C	16	76	62	20 (6.7)/7.7
G32581	ZFEN-D-16-62-91-C	16	91	62	20 (6.7)/7.7
G32585	ZFEN-D-16-62-106-C	16	106	62	20 (6.7)/7.7
G32589	ZFEN-D-16-62-121-C	16	121	62	20 (6.7)/7.7
G32553	ZFEN-D-20-28-76-C	20	76	28	20 (6.7)/7.7
G32556	ZFEN-D-20-28-91-C	20	91	28	20 (6.7)/7.7
G32558	ZFEN-D-20-28-106-C	20	106	28	20 (6.7)/7.7
G32560	ZFEN-D-20-28-121-C	20	121	28	20 (6.7)/7.7
G32563	ZFEN-D-20-45-76-C	20	76	45	20 (6.7)/7.7
G32567	ZFEN-D-20-45-91-C	20	91	45	20 (6.7)/7.7
G32571	ZFEN-D-20-45-106-C	20	106	45	20 (6.7)/7.7
G32575	ZFEN-D-20-45-121-C	20	121	45	20 (6.7)/7.7
G32578	ZFEN-D-20-62-76-C	20	76	62	20 (6.7)/7.7
G32582	ZFEN-D-20-62-91-C	20	91	62	20 (6.7)/7.7
G32586	ZFEN-D-20-62-106-C	20	106	62	20 (6.7)/7.7

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.


Distal Bifurcated Body Grafts (Continued)


Order Number	Reference Part Number	Ipsilateral Leg Diameter mm	Distal Body Length mm	Distal Leg Length mm	Introducer Sheath ID/OD Fr (mm)/mm
Distal Bifurcated Body Graft					
G32590	ZFEN-D-20-62-121-C	20	121	62	20 (6.7)/7.7
G32564	ZFEN-D-24-45-76-C	24	76	45	20 (6.7)/7.7
G32568	ZFEN-D-24-45-91-C	24	91	45	20 (6.7)/7.7
G32572	ZFEN-D-24-45-106-C	24	106	45	20 (6.7)/7.7
G32576	ZFEN-D-24-45-121-C	24	121	45	20 (6.7)/7.7
G32579	ZFEN-D-24-62-76-C	24	76	62	20 (6.7)/7.7
G32583	ZFEN-D-24-62-91-C	24	91	62	20 (6.7)/7.7
G32587	ZFEN-D-24-62-106-C	24	106	62	20 (6.7)/7.7
G32591	ZFEN-D-24-62-121-C	24	121	62	20 (6.7)/7.7

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

H&L-B One-Shot® Introduction System for Distal Body Grafts

20 Fr (6.7 mm) ID/7.7 mm OD






Maintains patency
and conformability
in a variety of iliac
anatomy.

Zenith[®] Spiral-Z[®]

AAA ILIAC LEG GRAFT

Zenith Flex
Zenith Renu
Zenith Fenestrated

Indications for Use



Indications for Use	For use with the Zenith® AAA Endovascular Graft family of products, during either a primary or a secondary procedure in patients who have adequate iliac/femoral access compatible with the required introduction systems.		
Ipsilateral Iliac Fixation Site	Distal fixation site >10 mm in length and 7.5-20 mm in diameter (measured outer wall to outer wall)		
Contralateral Iliac Fixation Site	Distal fixation site >10 mm in length and 7.5-20 mm in diameter (measured outer wall to outer wall)		
Access	Adequate femoral/iliac access compatible with the required introduction systems		
Introducer Diameters	9-16 mm grafts: ID 14 Fr (4.7 mm) OD 5.4 mm	20-24 mm grafts: ID 16 Fr (5.3 mm) OD 6.0 mm	
Contraindications	<ul style="list-style-type: none"> • Patients with known sensitivities or allergies to stainless steel, nitinol, polyester, solder (tin, silver), polypropylene, or gold • Patients with a systemic infection who may be at increased risk of endovascular graft infection 		

Contralateral/Ipsilateral Iliac Leg Grafts

Order Number	Reference Part Number	Iliac Leg Distal Diameter mm	Iliac Leg Working Length mm	Introducer Sheath ID/OD Fr (mm)/mm
9 mm Diameter				
G55218	ZSLE-9-39-ZT	9	39	14 (4.7)/5.4
G55219	ZSLE-9-56-ZT	9	56	14 (4.7)/5.4
G55220	ZSLE-9-74-ZT	9	74	14 (4.7)/5.4
G55221	ZSLE-9-90-ZT	9	90	14 (4.7)/5.4
G55222	ZSLE-9-107-ZT	9	107	14 (4.7)/5.4
G55223	ZSLE-9-122-ZT	9	122	14 (4.7)/5.4
11 mm Diameter				
G55224	ZSLE-11-39-ZT	11	39	14 (4.7)/5.4
G55225	ZSLE-11-56-ZT	11	56	14 (4.7)/5.4
G55226	ZSLE-11-74-ZT	11	74	14 (4.7)/5.4
G55227	ZSLE-11-90-ZT	11	90	14 (4.7)/5.4
G55228	ZSLE-11-107-ZT	11	107	14 (4.7)/5.4
G55229	ZSLE-11-122-ZT	11	122	14 (4.7)/5.4
13 mm Diameter				
G55230	ZSLE-13-39-ZT	13	39	14 (4.7)/5.4
G55231	ZSLE-13-56-ZT	13	56	14 (4.7)/5.4
G55232	ZSLE-13-74-ZT	13	74	14 (4.7)/5.4
G55233	ZSLE-13-90-ZT	13	90	14 (4.7)/5.4
G55234	ZSLE-13-107-ZT	13	107	14 (4.7)/5.4
G55235	ZSLE-13-122-ZT	13	122	14 (4.7)/5.4
16 mm Diameter				
G55236	ZSLE-16-39-ZT	16	39	14 (4.7)/5.4
G55237	ZSLE-16-56-ZT	16	56	14 (4.7)/5.4
G55238	ZSLE-16-74-ZT	16	74	14 (4.7)/5.4
G55239	ZSLE-16-90-ZT	16	90	14 (4.7)/5.4
20 mm Diameter				
G55240	ZSLE-20-39-ZT	20	39	16 (5.3)/6.0
G55241	ZSLE-20-56-ZT	20	56	16 (5.3)/6.0
G55242	ZSLE-20-74-ZT	20	74	16 (5.3)/6.0
G55243	ZSLE-20-90-ZT	20	90	16 (5.3)/6.0

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Contralateral/Ipsilateral Iliac Leg Grafts (Continued)


Order Number	Reference Part Number	Iliac Leg Distal Diameter mm	Iliac Leg Working Length mm	Introducer Sheath ID/OD Fr (mm)/mm
24 mm Diameter				
G55244	ZSLE-24-39-ZT	24	39	16 (5.3)/6.0
G55245	ZSLE-24-56-ZT	24	56	16 (5.3)/6.0
G55246	ZSLE-24-74-ZT	24	74	16 (5.3)/6.0
G55247	ZSLE-24-90-ZT	24	90	16 (5.3)/6.0

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Z-Trak[®] Introduction System for Iliac Legs

9-16 mm grafts: 14 Fr (4.7 mm) ID/5.4 mm OD

20-24 mm grafts: 16 Fr (5.3 mm) ID/6.0 mm OD



Distal Extensions

Order Number	Reference Part Number	Aortic Vessel OD mm	Diameter mm	Length mm	Introducer Sheath ID/OD Fr (mm)/mm
22 mm Diameter					
G35506	ZTA-DE-22-104-W	18-19	22	104	16 (5.3)/6.0
26 mm Diameter					
G35511	ZTA-DE-26-104-W	22-23	26	104	16 (5.3)/6.0
30 mm Diameter					
G35515	ZTA-DE-30-108-W	26-27	30	108	16 (5.3)/6.0
34 mm Diameter					
G35519	ZTA-DE-34-112-W	30	34	112	18 (6.0)/7.1
38 mm Diameter					
G35523	ZTA-DE-38-91-W	33-34	38	91	18 (6.0)/7.1
42 mm Diameter					
G36346	ZTA-DE-42-94-W	37-38	42	94	20 (6.7)/7.7
46 mm Diameter					
G38027	ZTA-DE-46-97-W	40-42	46	97	20 (6.7)/7.7

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Introduction System for Distal Extensions

22-30 mm grafts: 16 Fr (5.3 mm) ID/6.0 mm OD

34-38 mm grafts: 18 Fr (6.0 mm) ID/7.1 mm OD

42-46 mm grafts: 20 Fr (6.7 mm) ID/7.7 mm OD

